
Vasofix® Safety

B | BRAUN

B. Braun Melsungen AG
34209 Melsungen
Germany
www.bbraun.com

[LLDorder 2764 - G 150716](#)

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






Materials used

PP, PE, ABS, silicone rubber, chrome-nickel steel.

Catheter: FEP or Polyurethan see package labeling.

Indications

Creation of a secure peripheral venous access with a self-activating anti-needlestick device.

Blood transfusions or infusion of I. V. solutions suitable for administration via peripheral veins.

Intermittent intravenous drug administration.

Prophylactic creation of a secure venous access in patients who may require urgent intravenous drug administration, in particular prior to diagnostic or therapeutic procedures.

Contraindications

Vasofix® Safety should not be used in patients with known hypersensitivity to any of the materials employed.

Risks

This I. V. catheter is designed to reduce the risk of accidental needlesticks; however, care must be taken to avoid needlesticks. Universal precautions must be adhered to in accordance with Centers for Disease Control and Prevention/Occupational Safety and Health Administration (CDC/OSHA) standards for bloodborne pathogens, when starting or maintaining any I. V. catheter, to avoid the risk of exposure to contaminated blood.

Depending on how long the cannula is left in situ, on the type and amount of infusions or injections administered,

and on individual predisposition, thrombophlebitis may occur in the accessed vein.

Duration of use

Change according to CDC Guidelines and / or Hospital or Institutional protocols. The puncture site should be checked at regular intervals. Vasofix® Safety should be removed in the event of local or systemic signs of infection.

Warning

Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

After withdrawal, do not reintroduce the steel needle into the catheter, as the lat-

For single use only

See instructions for use

LOT

Batch number

ter may be cut off, leading to catheter embolism. Use only if packaging is intact.

Application

- 1 After disinfection of the puncture site and removal of the protective cap, puncture a suitable vein. If venipuncture was successful, blood will immediately be visible inside the transparent grip part.
- 2 Advance the catheter further into the vein, while slightly withdrawing the steel needle.
- 3 Using adhesive tape, fix catheter to the skin. The steel needle still in situ minimizes spillage of blood.
- 4 Before removing the steel needle compress the vein at the tip of catheter with the middle finger, to prevent spillage of blood. At the same time stabilize the catheter hub with the index finger to prevent catheter dislodgement during needle removal. Remove needle by pulling needle straight back. Metal safety clip will automatically attach to needle tip as needle tip exits catheter hub. Dispose of needle immediately into sharps container.
- 5 Connect to infusion line and cover puncture site with a sterile dressing.
- 6 Drugs can be injected without a cannula via the integrated injection port. To do this attach syringe conus with a slight twisting movement. Port opens and closes automatically.
- 7 Close protective cap of injection port.
- 8 For long-term interruption of I. V. administration introduce sterile mandrin (stylet). To avoid risk of any blood escaping press vein just beyond the tip of the catheter.

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